

ALTERNATE SEQUENCE POLISHING PADS



- Most polishing pads come in the sequence 50, 100, 200, 400, 800, 1500, 3000. The TST Alternate Sequence Pads run and perform differently. They are in 80, 150, 300, 500, 1000, 2000, 3000, 8500. This gives you two extra pads in the polishing part of the sequence.
- Many materials, like Blue Pearl and Black Galaxy, require a higher level of polishing to match the factory surface. These pads provide the extra polishing needs of specialty material surfaces.



GRANITE

QUARTZITE

Part Number	Size/Description
TSTASP-80	4", 80, Blue backing
TSTASP-150	4", 150, Yellow backing
TSTASP-300	4", 300, Orange backing
TSTASP-500	4", 500, Red backing
TSTASP-1000	4", 1000, Green backing
TSTASP-2000	4", 2000, Light blue backing
TSTASP-3000	4", 3000, Brown backing
TSTASP-8500	4", 8500, Black backing
TSTAPP-SET	Full Set